



PRECISA
INTERNATIONAL

Invitro Diagnostics

COMPACT ELISA KITS

www.PrecisaInternational.com | 305.669.0898



COMPACT ELISA



Cat. No.	Product Description	Total incubation time	Sensitivity	Packing	Notes
INFECTIOUS DISEASES					
E-101196	HIV-1/2 O	60 min (30/20/10)	>99.9%	96T	Qual.
E-128196	HIV-1/2 O Antigen/Antibody	120 min (60/30/30)	>99.9%	96T	Qual.
E-129196	HTLV I/II	90 min (30/30/30)	>99.9%	96T	Qual.
E-102196	HBsAg	70 min (60/10)	0.2ng/mL	96T	Qual.
E-103196	HCV	70 min (30/30/10)	>99.9%	96T	Qual.
E-104196	Syphilis (TP) IgG,IgM,IgA	105min (60/15/30)	99.9%	96T	Qual.
E-105196	HBsAb	45 min (30/15)	99.3%	96T	Qual.
E-127196	HbsAb Quantitative	45 min (30/15)	99.3%	96T	Quant.
E-106196	HBeAg	45 min (30/15)	99.7%	96T	Qual.
E-107196	HBeAb	45 min (30/15)	99.1%	96T	Qual.
E-108196	HBcAb TA (IgG+IgM+IgA)	45 min (30/15)	99.1%	96T	Qual.
E-121196	HEV IgM	75 min (30/30/15)	99.2%	96T	Qual.
E-123196	H pylori Antigen (Ag) (incl.extraction)	70 min (60/10)		96T	Qual./Quant
E-124196	H pylori IgG	70 min (30/30/10)	97.9%	96T	Qual./Quant
TORCH					
E-109196	Toxo IgG	70 min (30/30/10)	>99.9%	96T	Qual./Quant
E-110196	Toxo IgM	70 min (30/30/10)	>99.9%	96T	Qual.
E-111196	Rubella IgG	70 min (30/30/10)	96.4%	96T	Qual./Quant
E-112196	Rubella IgM	70 min (30/30/10)	90.0%	96T	Qual.
E-113196	CMV IgG	70 min (30/30/10)	>99.9%	96T	Qual./Quant
E-114196	CMV IgM	70 min (30/30/10)	89.3%	96T	Qual.
E-115196	HSV 1 IgG	70 min (30/30/10)	98.5%	96T	Qual.
E-116196	HSV1 IgM	70 min (30/30/10)	96.3%	96T	Qual
E-117196	HSV 2 IgG	70 min (30/30/10)	92.1%	96T	Qual.
E-118196	HSV 2 IgM	70 min (30/30/10)	90.5%	96T	Qual.
E-119196	HSV 1/2 IgG	70 min (30/30/10)	98.6%	96T	Qual.
E-120196	HSV 1/2 IgM	70 min (30/30/10)	85.7%	96T	Qual.
THYROID					
E-301196	TSH	75 min (60/15)	0.078uIU/mL	96T	Quant.
E-304196	Total T3	75 min (60/15)	0.1 ng/mL	96T	Quant.
E-302196	Total T4	75 min (60/15)	0.4 ug/dL	96T	Quant.
E-305196	Free T3	75 min (60/15)	0.24 pg/mL	96T	Quant.
E-303196	Free T4	75 min (60/15)	0.06 ng/dL	96T	Quant.
FERTILITY					
E-401196	Rapid-hCG	15 min (10/5)	0.6 mIU/mL	96T	Semi-Quant.
E-402196	LH	75 min (60/15)	0.5 mIU/mL	96T	Quant.
E-404196	Prolactin	75 min (60/15)	1.0 ng/mL	96T	Quant.
E-403196	FSH	75 min (60/15)	0.5 mIU/mL	96T	Quant.
TUMOR					
E-201196	AFP (Alpha-Fetoprotein)	45 min (30/15)	0.6 ng/mL	96T	Quant.
E-202196	CEA	45 min (30/15)	0.38 ng/mL	96T	Quant.
E-203196	PSA Total	45 min (30/15)	0.1 ng/mL	96T	Quant.

QUALITY PRODUCTS & SERVICES SINCE 1982

We offer up to date practical solutions for an ever changing medical device industry, with over 36 years of experience in the global markets. We inform. We link. We share in your success. Our goal is to share information on all products used in small to mid-size laboratories. We leverage our buying power and market knowledge to provide affordable, competitive goods to end users as well as distributors.

**TO ORDER PRODUCTS
CALL 305.669.0898
OR EMAIL
sales@PrecisaInternational.com**

