



**PRECISA**  
INTERNATIONAL

# Medical Equipment

## Microscopes

[www.PrecisaInternational.com](http://www.PrecisaInternational.com) | 305.669.0898

# MICROSCOPE COMPARISON CHART

## BEST SELLERS

## EPI FLUORESCENCE



	Innovation with BioVIEW	Innovation	i4	Revelation III	Lumin™ Innovation	Lumin™ i4
MODEL	Innovation with BioVIEW	Innovation	i4	Revelation III	Lumin™ Innovation	Lumin™ i4
MODEL NUMBER	INS-T4BV-IPL3*	INM-B04A-iPL3*	I4M-BN4A-IPL3	R3M-BN4A-DAL3	INS-EPB4-IPL3*	I4S-EPB4-IPL3
HEAD	Trinoc	Binoc (Trinoc available)	Binoc (Trinoc available)	Binoc (Trinoc available)	Binoc (Trinoc available)	Binoc (Trinoc available)
STAGE	Rackless Mechanical	Rackless Mechanical	Mechanical	Mechanical	Rackless Mechanical	Mechanical
OBJECTIVES	4-10-40-100 High Contrast Infinity Plan	4-10-40-100 High Contrast Infinity Plan	4-10-40-100 Infinity Plan	4-10-40-100 DIN Achromat	4-10-40-100 High Contrast Infinity Plan	4-10-40-100 Infinity Plan
NOSEPIECE	Reverse Quadruple	Reverse Quadruple	Reverse Quadruple	Quadruple Forward-Facing	Reverse Quadruple	Reverse Quadruple
FEATURES	HD Camera & HD Monitor with Trinoc InnoVation One Part Number - All Included	High Resolution, Flat-field, Rackless Stage, LED Light	High Resolution, 100% Flat-field, LED Light	Affordable, Durable, LED Light	Images rival those obtained from traditional, expensive epifluorescence systems	Images rival those obtained from traditional, expensive epifluorescence systems
APPLICATIONS	Research, Teaching, Demonstration, Pathology, Cytology, Oncology, Hematology	Pathology, Hematology, Oncology, Gynecology, Research, Microbiology	Cytology, Hematology , Oncology, Gynecology, Research, Microbiology	Routine Lab, POL, Tech Schools, Universities	Tuberculosis, Malaria, STD, Pathology, Immunology, Research	Tuberculosis, Malaria, STD, Pathology, Immunology, Research

\*5 Objective Model  
with 20X Available

\*5 Objective Model  
with 20X Available

\*5 Objective Model  
with 20X Available

# MICROSCOPE COMPARISON CHART

## SPECIALIZED



	i4 Semen Evaluation	Mi5 Polarizer	Mi5 Dual Head	Inverted
MODEL				
MODEL NUMBER	i4S-SEB4-iPL3	M5M-BPOL-IPH3	M5M-DT5A-IPL3	L1M-T04A-DPH3
HEAD	Binoc (Trinoc available)	Binoc (Trinoc available)	Binoc + Trinoc*	Trinoc
STAGE	Mechanical - Heated	Round, Rotates 360, Mechanical Optional	Mechanical	Mechanical
OBJECTIVES	4-10-40-100 DIN Phase Plan	4-10-20-40 Infinity Plan	4-10-20-40-100 Infinity Plan	4-40 Plan & 10-20 Plan Phase
NOSEPIECE	Reverse Quadruple	Reverse Quintuple	Reverse Quintuple	Reverse Quintuple
FEATURES	Heated Stage, Phase For Semen, Brightfield	High-resolution Optics, Polarizing Analyzer with 360 Degree Increments	Two People Viewing Same Sample	Centering Eyepiece for Phase Objective, Coaxial Drive Controls
APPLICATIONS	Fertility Clinics, IVF, Research	Rheumatology, Urology	Pathology, Hematology, Oncology, Gynecology, Research, Microbiology	Pathology, Biology, Medical Researcher

\*Binoc + Binoc Model Available

# QUALITY PRODUCTS & SERVICES SINCE 1982

We offer up to date practical solutions for an ever changing medical device industry, with over 36 years of experience in the global markets. We inform. We link. We share in your success. Our goal is to share information on all products used in small to mid-size laboratories. We leverage our buying power and market knowledge to provide affordable, competitive goods to end users as well as distributors.

**TO ORDER PRODUCTS  
CALL 305.669.0898  
OR EMAIL  
[sales@PrecisaInternational.com](mailto:sales@PrecisaInternational.com)**

